

Vibratory Polishing of High-Purity Metals to Improve EBSD Results

EBSD experiments evaluated the band contrast of eighteen (18) high-purity (generally >99.95%) specimens mechanically prepared using five steps (four for Ti). These specimens varied from Mg (atomic number 12) to Bi (atomic number 83) and covered the range of metallic crystal structures: body-centered cubic (6), face-centered cubic (4), hexagonal close-packed (5), diamond cubic (1) and rhombohedral/trigonal (2). The table below lists the specimens prepared using our standard methods and analyzed. Results for six of these after vibratory polishing are shown in the second table.


Specimens of pure Sb, V and Zr were susceptible to SiC embedment, even though the grit size was coarse, e.g., 240- and 320-grit. Hence, grinding was repeated after coating the paper with paraffin wax. Attack polishing was used, mainly with 30% conc. H₂O₂, for the last step for preparing Cr, Nb, Ti, W and Zr. MasterMet® colloidal silica was used for the last step, except for preparing Fe (MasterPrep® alumina was used) and Mg (water-free MasterPolish® was used). Oil-based diamond suspensions (9-, 3- and 1-µm) were used to prepare the high-purity (99.999%) Mg. For the Bi and Pb pure specimens, grinding used four steps: 240-, 320-, 400- and 600-grit SiC paper coated with paraffin wax with low loads, followed by three polishing steps using 5-, 1- and 0.3-µm alumina slurries and a final polish with MasterMet® colloidal silica. All polishing steps used MicroCloth® synthetic suede cloth. Although the Bi produced an excellent EBSD pattern, none was obtained with the pure Pb specimen. A one-hour vibratory polish with MasterMet® colloidal silica using a MicroCloth® pad was required to obtain a diffraction pattern for Pb.

The result for pure Zr in the table below was obtained on the same specimen as illustrated here, but after an earlier preparation attempt with a less effective preparation method than used to generate the grain maps. The average band contrast for the high-purity Zr specimen using the improved preparation method was 92.34 and ~90% of the pixels produced indexable diffraction patterns. For the initial effort with high-purity Zr, which yielded an average band contrast of 77.3, only about 20% of the pixels in the SEM field yielded indexable pixels. Preparation methods are improved when they do not perform adequately. Zr is among the most difficult metals to prepare for EBSD.


Two maps of high-purity (99.99%), annealed Zr are shown. The first was constructed by adding an all Euler grain map with a band contrast map; the second shows an inverse pole figure map figure map, plus grain boundaries, with the grains with missing pixels (black spots in the first map) filled in. The band contrast averaged 92.34 for the area shown.

Band Contrast Values for 18 Pure Metals

High-Purity Elements	Atomic Number	Crystal Structure	Band Contrast (0-255)
Mg	12	hcp	161.2
Al	13	fcc	151.2
Si	14	diamond cubic	205.75
Ti	22	hcp	134.0
V	23	bcc	102.2
Cr	24	bcc	88.27
Fe	26	bcc	105.4
Ni	28	fcc	85.0
Cu	29	fcc	122.6
Zn	30	hcp	170.8
Zr	40	hcp	77.3
Nb	41	bcc	145.6
Ru	44	hcp	66.0
Sb	51	rhombohedral	180.2
Ta	73	bcc	122.8
W	74	bcc	91.6
Pb	82	fcc	no pattern
Bi	83	rhomb./trigonal	255



All Euler Angles plus Band Contrast map for high-purity (99.99%) Zr




Inverse pole figure map with grain boundaries in black for high-purity (99.99%) Zr


Improved Preparation Method for Zirconium and Alloys

Surface	Abrasive	Lubricant	Force (lbs)	Speed (rpm)	Direction	Time (minutes)
SiC Paper (wax coated)	240-grit	Water	5	240	Contra	Until Planar
SiC Paper (wax coated)	320-grit	Water	5	240	Contra	1
UltraPol silk cloth	9-µm Diamond Paste	MetaDi Fluid	5	200	Contra	10
TriDent polyester cloth	3-µm Diamond Paste	MetaDi Fluid	5	200	Contra	7
TriDent polyester cloth	1-µm Diamond Paste	MetaDi Fluid	5	200	Contra	5
MicroCloth pad	0.05-µm Colloidal Silica*	None	6	200	Contra	7

* Attack polish with 1 part H₂O₂ (30%) to 5 parts Colloidal Silica



EBSD Pattern for annealed 99.99% Zr



EBSD Pattern for EB Melted (>99%) Zr (hot extruded and cold drawn)

Band Contrast Improvement Due to Vibratory Polishing (20 min.*)

High-Purity Element	Mean Band Contrast (0 to 255)	
	Standard Method	Standard + Vibratory Polish
Mg	161.2	175.25 (+8.7%)
Si (single crystal)	205.75	233 (+13.2%)
Ti	134.0	146.2 (+9.1%)
Ni	85.0	102.8 (+20.9%)
Nb	145.6	151.2 (+3.8%)
Pb	No pattern	108.0

* A 60 minute vibratory polish was used for the lead specimen.



Need Help Preparing EBSD Specimens? Call your Buehler representative for excellence in equipment, consumables, methods and training – for any crystalline material. See www.buehler.com and join our e-club for technical support.